

DATA SHEET FRONTLINE - Bleed Control Trauma Kit

SKU Number **SS-6010125040**
 EAN13 Barcode **5060660858163**

Date of Publication **09-01-2025**
 Commodity Code **30065000.00**



<p>Product Description</p>	<p>FRONTLINE – Bleed Control Trauma Kit</p> <p>When severe haemorrhage control is needed, having a kit equipped with the right medical products is essential. Simple, user-friendly Bleed Control Kits are critical to the effectiveness of first responders.</p> <p>The STEIN FRONTLINE Bleed Control Trauma Kit includes high-quality medical supplies necessary to manage life-threatening haemorrhages.</p> <p>This kit allows trained First Aiders, Advanced Professionals, and First Responders to provide lifesaving care in critical situations whilst waiting for Emergency Services to arrive.</p> <p>Contained in a durable, lightweight bag made from advanced materials, this kit is designed for quick access in emergencies.</p> <p>The STEIN FRONTLINE Bleed Control Trauma Kit functions as a grab bag that can be easily attached to the harness of a climber performing an Aerial Rescue. For added convenience, it can also be clipped near the casualty, allowing rapid access to its contents for immediate treatment.</p>
----------------------------	---

Contents	<ul style="list-style-type: none"> 1 x SAM XT Tourniquet 1 x 6" FCP Military Emergency Bandage 1 x ChitoSAM Haemostatic Gauze 4' Z-Fold 1 x Large First Aid Dressing 18cm x 18cm 3 x 20ml Ampoules Sterile Wound Wash 1 x High Performance Medical Shears 1 x Medical Marker Pen 2 x Disposable Medical Gloves 3 x Security Seals
-----------------	--

Safety Standards	<ul style="list-style-type: none"> • All medical products used in the STEIN Kits are CE Marked and conform with the relevant European Union harmonisation legislation.
-------------------------	---

Product Weight	<ul style="list-style-type: none"> • 679g
Size	<ul style="list-style-type: none"> • 270cm x 170cm x 70cm (approx.)

Storage	<ul style="list-style-type: none"> • Keep stored and transported at ambient temperature (0°C to + 30°C) • Do not store in direct/high heat or sunlight as this may affect the internal components of the kit
Lifetime	<ul style="list-style-type: none"> • Not all products have an expiry date, however most sterile items do. If a non-sterile product has no expiry date and still sealed in its original packaging, it's a question of judgement as to whether it is still suitable for use. • Items that have exceeded their expiry date should be disposed and replaced as they may not be effective and can start to deteriorate, making them unreliable.
Cleaning	<ul style="list-style-type: none"> • Clean the storage bag using warm soapy water and soft cloth
Disposal	<ul style="list-style-type: none"> • Recycling is preferred where possible. Look out for the recycling symbol for material categories

Replacement Components	<ul style="list-style-type: none">• Following the use of any product the kit should be replenished with items used so they are fully functional for future incidents• Refer to the STEIN Website for replacement items www.steinworldwide.com
Additional Information	<ul style="list-style-type: none">• In extreme cases of haemorrhage further products may be required over and above the items in these kits.• This kit does not contain products for the treatment of minor injuries• This kit does not replace the mandatory Workplace First Aid Kits complying with BS 8599-1:2019
Training	<ul style="list-style-type: none">• It is not feasible to address every potential scenario involving the use of this product. Users are strongly advised to obtain professional training from a fully qualified and competent instructor before use.• Training Providers can be recommended upon request
Trademarks	<ul style="list-style-type: none">• STEIN® is a registered trademark