

DATA SHEET SOLO - Personal First Aid Pack

SKU Number **SS-6032650030**
 EAN13 Barcode **5060483285337**

Date of Publication **09-01-2025**
 Commodity Code **30065000.00**

Product Image			
---------------	---	--	---

Product Description	<p>SOLO - Personal First Aid Pack</p> <p>The STEIN SOLO Personal First Aid Pack is an essential piece of equipment for every arborist. Compact and practical, it equips you to handle a wide range of minor injuries commonly encountered in the arboricultural industry.</p> <p>Packed with necessary essentials, this kit ensures you're prepared to treat small lacerations and minor injuries.</p> <p>Contained in a robust, weather-resistant pouch which features a sealed zipper for added protection. It is equipped with a Velcro belt loop with a securing press stud, and a D-ring, allowing it to attach securely to any climbing harness for convenient access.</p>
---------------------	--

Contents	<ul style="list-style-type: none"> 1 x Medium Wound Dressing No. 14 1 x Eye/Wound Solution (20ml Ampoule) 1 x Pack of Plasters (10) 1 x Pack of Sterile Skin Closure Strips (5) 1 x Disposable Medical Gloves 2 x Sterile Wound Cleansing Wipes 1 x First Aid Guidance Leaflet
-----------------	---

Safety Standards	<ul style="list-style-type: none"> • All medical products used in the STEIN Kits are CE Marked and conform with the relevant European Union harmonisation legislation.
-------------------------	---

Product Weight	<ul style="list-style-type: none"> • 153g
Size	<ul style="list-style-type: none"> • 100cm x 140cm x 50cm (approx.)

Storage	<ul style="list-style-type: none"> • Keep stored and transported at ambient temperature (0°C to + 30°C) • Do not store in direct/high heat or sunlight as this may affect the internal components of the kit
Lifetime	<ul style="list-style-type: none"> • Not all products have an expiry date, however most sterile items do. If a non-sterile product has no expiry date and still sealed in its original packaging, it's a question of judgement as to whether it is still suitable for use. • Items that have exceeded their expiry date should be disposed and replaced as they may not be effective and can start to deteriorate, making them unreliable.
Cleaning	<ul style="list-style-type: none"> • Clean the storage bag using warm soapy water and soft cloth
Disposal	<ul style="list-style-type: none"> • Recycling is preferred where possible. Look out for the recycling symbol for material categories

Replacement Components	<ul style="list-style-type: none">• Following the use of any product the kit should be replenished with items used so they are fully functional for future incidents• Refer to the STEIN Website for replacement items www.steinworldwide.com
Additional Information	<ul style="list-style-type: none">• In extreme cases of injury further products may be required over and above the items in these kits.• This kit is not suitable to treat a severe life-threatening haemorrhage• This kit does not replace the mandatory Workplace First Aid Kits complying with BS 8599-1:2019
Training	<ul style="list-style-type: none">• It is not feasible to address every potential scenario involving the use of this product. Users are strongly advised to obtain professional training from a fully qualified and competent instructor before use.• Training Providers can be recommended upon request
Trademarks	<ul style="list-style-type: none">• STEIN® is a registered trademark